## 新型恆溫式圈環形核酸增幅法簡介與應用 Loop-mediated Isothermal Amplification (LAMP)

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- 1、Notomi T, Okayama H, Masubuchi H et al. Loop-mediated isothermal amplification of DNA. *Nucleic Acids Res.* 2000;28:E63.
- 2、Nagamine K, Hase T, Notomi T. Accelerated reaction by loop-mediated isothermal amplification using loop primers. *Mol. Cell Probes* 2002;16:223-229.
- 3、Iwamoto T, Sonobe T, Hayashi K. Loop-mediated isothermal amplification for direct detection of *Mycobacterium tuberculosis* complex, *M. avium*, and *M. intracellulare* in sputum samples. *J. Clin. Microbiol.* 2003;41:2616-2622.
- 4、Mori Y, Nagamine K, Tomita N et al. Detection of loop-mediated isothermal amplification reaction by turbidity derived from magnesium pyrophosphate formation. *Biochem. Biophys. Res. Commun.* 2001;289:150-154.
- 5、Maruyama F, Kenzaka T, Yamaguchi N et al. Detection of bacteria carrying the *stx2* gene by in situ loop-mediated isothermal amplification. *Appl. Environ. Microbiol.* 2003;69:5023-5028.
- 6、Fukuta S, Iida T, Mizukami Y et al. Detection of Japanese yam mosaic virus by RT-LAMP. *Arch. Virol.* 2003;148:1713-1720.
- 7、Parida M, Posadas G, Inoue S et al. Real-time reverse transcription loop-mediated isothermal amplification for rapid detection of West Nile virus. *J. Clin. Microbiol.* 2004;42:257-263.
- 8、Ihira M, Yoshikawa T, Enomoto Y et al. Rapid diagnosis of human herpesvirus 6 infection by a novel DNA amplification method, loop-mediated isothermal amplification. *J. Clin. Microbiol.* 2004;42:140-145.
- 9、Kuboki N, Inoue N, Sakurai T et al. Loop-mediated isothermal amplification for detection of African trypanosomes. *J. Clin. Microbiol.* 2003;41:5517-5524.
- 10、Fukuta S, Kato S, Yoshida K et al. Detection of tomato yellow leaf curl virus by loop-mediated isothermal amplification reaction. *J. Virol. Methods* 2003;112:35-40.
- 11、Kono T, Savan R, Sakai M et al. Detection of white spot syndrome virus in shrimp by loop-mediated isothermal amplification. *J. Virol. Methods* 2004;115:59-65.

- 12、Savan R, Igarashi A, Matsuoka S et al. Sensitive and rapid detection of edwardsiellosis in fish by a loop-mediated isothermal amplification method. *Appl. Environ. Microbiol.* 2004;70:621-624.
- 13、Nagamine K, Kuzuhara Y, Notomi T. Isolation of single-stranded DNA from loop-mediated isothermal amplification products. *Biochem. Biophys. Res. Commun.* 2002;290:1195-1198.
- 14、Poon LL, Leung CS, Tashiro M et al. Rapid Detection of the Severe Acute Respiratory Syndrome (SARS) Coronavirus by a Loop-Mediated Isothermal Amplification Assay. *Clin. Chem.* 2004.

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